P/N: 110521-U

## Battery type: NiMH

**Kinetic Concepts Inc**  
**VAC ATS Wound Vacuum**  
**Battery Typ M3252795-A**  
**12 Volt 3,7 Ah**

### User Instructions

#### **Installation and use**

This product must be installed by a qualified professional!  
Refer to the equipment manufacturer's instructions!"

Before its first use, fully charge the battery using the specified charger. The battery may only be replaced by authorised technical personnel with the device switched off. Avoid mechanical impact on the battery and its terminals. Do not allow the terminals to come into contact with metallic objects:

#### **RISK OF SHORT CIRCUIT, FIRE AND EXPLOSION**

Do not expose the battery to fire or high temperatures. Do not open the battery.

#### **Charging the battery**

Before charging the battery, ensure correct polarity! To operate at its specified capacity, always charge the battery in room temperature. The battery may heat up while charging. Charge the battery at the rate and for the time specified. Avoid total discharge of the battery, i.e. do not leave the battery for longer periods in the unused, activated device. Always make sure to fully charge NiMH (Nickel, Metal Hydrid) batteries to avoid loss of capacity through the lazy battery effect. Some batteries only achieve their full capacity after several charging cycles. We therefore recommend charging the battery in the device before its first use and then discharging and recharging it. This will increase the useful life of the battery.

#### **Storing, maintaining the battery**

Store the battery in a cool, dry and well-ventilated place, protected from damp and humidity. Protect from short circuits. If several products are included in the pack, leave the instructions in the outer packaging until all products are used.

The batteries are provided with a LOT no. or production code: 1234, where 1=day of week, 2+3=calendar week; 4=year.  
Example: 2455 = Tue, CW 45, 2015.

Stored batteries should be recharged at least every 3 months. Batteries already in use should be checked for efficiency once a year.

Batteries achieving less than 80% of their rated capacity should be discarded or replaced. Batteries older than 24 months should be checked every 6 months.

#### **Disposal of used batteries**

Used batteries must not be disposed of with common household waste! Please take old batteries to a designated collection point. Details and information on disposal are available at:  
<http://www.grs-batterien.de>

Should you have any queries, please contact us at any time using the following free service hotline:

**0800-AKKUMED** (0800-2558633).

Akkuplanet GmbH; Berta-Benz-Str. 18; D - 40670 Meerbusch;  
Germany; phone +49-2159-922 40-00

CE mark based on the Medical Devices Directive 93/42/EEC of 14 June 1993

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